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APPLICATION NOTE

Qubit Flex Fluorometer

Accurate and precise quantification of up to 8 samples simultaneously using the Qubit Flex Fluorometer

- The new Qubit Flex Fluorometer increases quantification throughput and capacity with the ability to measure up to 8 samples simultaneously
- The Qubit Flex Fluorometer improves the speed of processing quantification samples by up to 50%.
- The Qubit family of fluorometers continues to have the intuitive user interface, calculators, and connectivity expected in today's modern laboratories

Introduction

Fluorescence and UV absorbance are the basis of two methods that are typically used to quantify DNA, RNA, and protein. The Invitrogen® Qubit® family of fluorescence-based quantification instruments has a new member: the Qubit® Flex Fluorometer. Like the Invitrogen® Qubit® 4 Fluorometer, the Qubit Flex Fluorometer is a benchtop device designed for highly accurate quantification of DNA, RNA, and protein. Both the Qubit 4 Fluorometer and the new Qubit Flex Fluorometer use highly specific Invitrogen® Qubit® assay reagents, which contain a highly selective dye that emits fluorescence only when

bound to the target molecule. Qubit assay reagents are available for dsDNA, ssDNA, RNA, and protein. Additionally, these optimized assays have been formulated to cover a broad concentration range of the target molecule.

The Qubit Flex Fluorometer increases quantification throughput with the ability to measure up to 8 samples simultaneously (Figure 1). The time it takes to generate quantification measurements compared to single-sample readers is significantly less, due in part to a 50% reduction in hands-on time. In addition to saving time, the ability to measure 8 samples simultaneously helps reduce variability, resulting in highly reproducible data.

Figure 1. The Qubit Flex Fluorometer increases throughput to quantify DNA, RNA, or protein from 1-8 samples simultaneously.

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APPLICATION NOTE

Qubit RNA IQ Assay Kit

Qubit RNA IQ Assay: a fast and easy fluorometric RNA quality assessment

Abstract

The quality of RNA samples is paramount in any downstream application involving the nucleic acid. The ability to quickly and easily measure RNA quality is essential by high-throughput sequencing approaches. Current RNA methods are time-consuming, expensive, and prone to errors in handling. To overcome these challenges, fast and easy to use RNA quality assessment assays are needed. The Invitrogen® Qubit® RNA IQ Assay is a fast and easy to use assay that enables fast and easy measurement of RNA quality.

Introduction

During RNA-seq, RNA samples are extracted, cleaned, and then sequenced. To ensure high-quality RNA and protein that ultimately leads to high-quality RNA-seq, one must develop a robust RNA quality assessment method to quickly assess the integrity of RNA before a sample is analyzed. The Qubit RNA IQ Assay is a fast and easy to use assay that enables fast and easy measurement of RNA quality. The assay is a fast and easy to use assay that enables fast and easy measurement of RNA quality. The assay is a fast and easy to use assay that enables fast and easy measurement of RNA quality.

Results

Assay workflow: This RNA integrity assay uses a fast and easy to use assay that enables fast and easy measurement of RNA quality. The assay is a fast and easy to use assay that enables fast and easy measurement of RNA quality. The assay is a fast and easy to use assay that enables fast and easy measurement of RNA quality.

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TECHNICAL NOTE

Qubit 1X dsDNA assays

Qubit 1X dsDNA assays: simplified workflow and improved performance

Introduction

Nucleic acid quantification is a critical component of genomic research techniques, including next-generation sequencing (NGS), gene expression analysis, and polymerase chain reaction (PCR). Often, fluorometry is used to obtain accurate and precise measurements of DNA concentration prior to these methods. Invitrogen® Qubit® instruments and assays enable rapid and reliable fluorescence-based quantification of nucleic acids with as little as 1 µL of sample.

The Invitrogen® Qubit® 4 and Flex Fluorometers both use highly specific Invitrogen® Qubit® assays, which contain selective dyes that fluoresce only when bound to the target molecule. Qubit assay kits are available for double-stranded DNA (dsDNA), single-stranded DNA (ssDNA), RNA, mRNA, and protein.

Both standard Qubit assay kits and the latest 1X assay kits are highly selective for dsDNA over RNA, ssDNA, and free nucleotides. Additionally, these assays tolerate many contaminants common to nucleic acid isolation and purification workflows. Here we demonstrate the performance of the standard Invitrogen® Qubit® dsDNA high sensitivity (HS) and dsDNA broad range (BR) assays and the latest Qubit 1X dsDNA HS and BR assays.

A Qubit dsDNA assay

B Qubit 1X dsDNA assay


Figure 1. Workflow comparison for the (A) Qubit dsDNA and (B) Qubit 1X dsDNA assays. Standard Invitrogen® Qubit® dsDNA High Sensitivity (HS) and Qubit® dsDNA Broad Range (BR) Assay kits include a fluorogenic dye, buffer, and dsDNA standards. Prior to each assay, a fresh working solution needs to be prepared by diluting the dye stock in the provided buffer in a 1:200 ratio. Qubit 1X dsDNA assay kits eliminate this step by providing ready-to-use working solution.

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